

# Efficiency in Application

All in one

BeneFusion eSP and eVP satisfy various infusion purposes by combining all functions together.



## Integrated central monitoring

BeneVision CMS Vet offers one-stop monitoring of all patients' vital sign and infusion treatment details, providing comprehensive information for veterinarian to improve the quality of patient care.



[www.mindray.com](http://www.mindray.com)

P/N:ENG-XXXX-210285X12P-20200805  
©2020 Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All rights reserved.

**mindray**  
healthcare within reach



**mindray**

## BeneFusion eSP/eVP/eDS Vet

Infusion System

# Efficiency in every droplet




# Efficiency in Workflow

## Fast preparation

### SmartRapid™

SmartRapid™ ensures timely infusion by significantly shortening the start-up time, from turning on the pump to delivering the first drop of medication.




## Intuitive interaction

- 3.5" Colored capacitive touchscreen brings users an excellent experience with smooth operation.
- Intuitive user interface enables quick programming of key parameters, making workflow process much easier.



## Flexible infusion adjustment


Easily adjusting the infusion therapy within 3 steps.



# Efficiency in Safety

## Precise infusion performance

- High accuracy: eSP  $\pm 1.8\%$ ; eVP  $\pm 5\%$
- Long-hour accuracy assured



## SmartAIR™

With SmartAIR™, BeneFusion eVP significantly enhances the safety level of IV administration.

- Dual ultrasound sensor to detect the air bubbles more precisely, avoiding missing or false air-in-line alarm.
- 15ul air bubble detection size on BeneFusion eVP, ensuring patient safety throughout the infusion, even for neonates.



## Color coding

- The color coding of drug name assists users to easily select and verify the correct drug.



## Quick problem solving

### Dynamic Pressure System

Speedometer style indicator with numerical pressure to monitor the in-line pressure trend just with a glance.



Alert with a visual message about the possible occlusion before interrupting the infusion.

Auto-attempting mechanism enables infusion resume as early as possible to ensure continuous infusion.

### Intuitive Alarm System

Instructional animation pops up to guide users to quickly solve the problem.

